

Ultra accuracy and safety in venous blood collection

Blood collection needle

Easy puncture, minimal pain

- I Also called multi-sampling needle, which is strictly designed and manufactured according to EC directive 93/42/EEC, and used with vacuum blood collection tube and plastic holder for venous blood collection.
- I Consist of cannula, rubber protector, top protector bottom protector and hub
- I Medical grade stainless steel cannula.
- I EO sterile, Non toxic, non pyrogenic
- I Country of Origin: China
- I Validity: 5 years

CE0197



| Cat No. | Product | Gauge x Length | Color | Packing info |
|---------|---------------------------------------|----------------|--------|-------------------------|
| N-N3M | Blood collection needle, sterile, 20G | 20G x 1 1/2" | Yellow | 100pcs/box, 5000pcs/ctn |
| N-N3P | Blood collection needle, sterile, 21G | 21G x 1 " | Green | 100pcs/box, 5000pcs/ctn |
| N-N3Q | Blood collection needle, sterile, 21G | 21G x 1 1/2" | Green | 100pcs/box, 5000pcs/ctn |
| N-N3S | Blood collection needle, sterile, 22G | 22G x 1" | Black | 100pcs/box, 5000pcs/ctn |
| N-N3T | Blood collection needle, sterile, 22G | 22G x 1 1/2" | Black | 100pcs/box, 5000pcs/ctn |

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

www.nubenco.com.cn