

Blood lancets

Easy puncture, minimal pain

- I Fine gauge, tri-bevel tip for virtually painless sampling
- I Universal design fits almost all lancing devices
- I Disposable, stainless steel type, with fine point for controlled penetration
- Blood lancet with twist cap, lancet devices available
- I Medical grade stainless steel needle.
- I Gamma Radiation sterilization, Non toxic, non pyrogenic
- I Country of Origin: China
- I Validity: 5 years
- I Packing: blister packed for stainless steel type

CE0197



Cat No.	Product	Gauge x Length	Packing info
N-N6	Blood Lancet with twist cap, sterile, blue	28G	200pcs/box, 20000pcs/ctn
N-N64	Blood Lancet with twist cap, sterile, light blue	28G x 1/2"	100pcs/box, 5000pcs/ctn
N-N7	Blood Lancet (Scarifier), sterile, disposable		200pcs/box, 20000pcs/ctn
N-N8	Lancet holder (pen-type), regular size		box x 1pcs, 100pcs/ctn

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

www.nubenco.com.cn