

Pasteur Pipettes (Transfer)

- I Disposable plastic pipette
- I Graduated 1ml, 3ml
- I Laboratory use, for liquid transfer
- I Non sterile
- I Country of Origin: China
- I Validity: 5 years



- Ø 1ml, 3ml graduated
- Ø Box of 500pcs, 5000pcs/carton

Cat. No.	Product	Graduation	Carton Size	Packing info
N-EB	Pasteur pipette, transfer	1ml	78mm x 44mm x 29mm	500pcs/box, 5000pcs/ctn
N-EC	Pasteur pipette, transfer	3ml	80mm x 59mm x 23mm	500pcs/envelop, 5000pcs/ctn

Mfg.: Nubenco (Shanghai) Co, Ltd
Shanghai, China
Tel.: +86-21-58967745
www.nubenco.com.cn