

Ultra accuracy in dosage and safety

Hypodermic needle

Easy puncture, minimal pain

CE 0197



- Ø Fits luer slip and luer lock connector
- Ø High grade stainless steel
- Ø Sharp bevel and siliconized cannula
- Ø Various gauges 16-30G, length 1/2", 1", 1 1/2" available
- Ø EO sterile, individual blister packed
- Ø Validity: 5 years

Cat No.	Product	Color	Gauge x Length	Packing
N-N1B	Hypodermic needle, sterile, blister packed	White	16G x 1/2"	Box x 100pcs, 10000pcs/carton
N-N1A	Hypodermic needle, sterile, blister packed	White	16G x 1"	Box x 100pcs, 10000pcs/carton
N-N16	Hypodermic needle, sterile, blister packed	White	16G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1F	Hypodermic needle, sterile, blister packed	Pink	18G x 1"	Box x 100pcs, 10000pcs/carton
N-N1G	Hypodermic needle, sterile, blister packed	Pink	18G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1J	Hypodermic needle, sterile, blister packed	Beige	19G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1I	Hypodermic needle, sterile, blister packed	Beige	19G x 1"	Box x 100pcs, 10000pcs/carton
N-N1L	Hypodermic needle, sterile, blister packed	Yellow	20G x 1"	Box x 100pcs, 10000pcs/carton

N-N1M	Hypodermic needle, sterile, blister packed	Yellow	20G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1Q	Hypodermic needle, sterile, blister packed	Green	21G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1P	Hypodermic needle, sterile, blister packed	Green	21G x 1"	Box x 100pcs, 10000pcs/carton
N-N1S	Hypodermic needle, sterile, blister packed	Black	22G x 1 "	Box x 100pcs, 10000pcs/carton
N-N1T	Hypodermic needle, sterile, blister packed	Black	22G x 1 1/2 "	Box x 100pcs, 10000pcs/carton
N-N1V	Hypodermic needle, sterile, blister packed	Blue	23G x 1"	Box x 100pcs, 10000pcs/carton
N-N1W	Hypodermic needle, sterile, blister packed	Blue	23G x 1 1/2"	Box x 100pcs, 10000pcs/carton
N-N1Y	Hypodermic needle, sterile, blister packed	Purple	24G x 1"	Box x 100pcs, 10000pcs/carton
N-N1AE	Hypodermic needle, sterile, blister packed	Orange	25G x 1"	Box x 100pcs, 10000pcs/carton
N-N11	Hypodermic needle, sterile, blister packed	Orange	25G x 5/8"	Box x 100pcs, 10000pcs/carton
N-N12	Hypodermic needle, sterile, blister packed	Brown	26G x 1/2"	Box x 100pcs, 10000pcs/carton
N-N13	Hypodermic needle, sterile, blister packed	Grey	27G x 1/2"	Box x 100pcs, 10000pcs/carton
N-N1AL	Hypodermic needle, sterile, blister packed	Salmon	29G x 1/2"	Box x 100pcs, 10000pcs/carton
N-N15	Hypodermic needle, sterile, blister packed	Ice	30G x 1/2"	Box x 100pcs, 5000pcs/carton

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

Fax: +86-21-58967742

Email:info@nubenco.com.cn

www.nubenco.com.cn