

**Ultra accuracy and safety in venous blood collection**  
**Easy puncture, minimal pain**

**Blood collection needle**

**Multi-sampling needle**

- I Strictly designed and manufactured according to EC directive 93/42/EEC, and used with vacuum blood collection tube and plastic holder for venous blood collection.
- I Consist of cannula, rubber protector, top protector bottom protector and hub
- I Medical grade stainless steel cannula.
- I EO sterile, Non toxic, non pyrogenic
- I Country of Origin: China
- I Validity: 3 years

**CE0197**



Cat No.	Product	Gauge x Length	Color	Packing info
N-N3M	Blood collection needle, multi-sampling type, sterile	20G x 1 1/2"	Yellow	100pcs/box, 5000pcs/ctn
N-N3P	Blood collection needle, multi-sampling type, sterile	21G x 1 "	Green	100pcs/box, 5000pcs/ctn
N-N3Q	Blood collection needle, multi-sampling type, sterile	21G x 1 1/2"	Green	100pcs/box, 5000pcs/ctn
N-N3S	Blood collection needle, multi-sampling type, sterile	22G x 1"	Black	100pcs/box, 5000pcs/ctn
N-N3T	Blood collection needle, multi-sampling type, sterile	22G x 1 1/2"	Black	100pcs/box, 5000pcs/ctn

**Ultra accuracy and safety in venous blood collection**  
**Easy puncture, minimal pain**

**Blood collection needle**

**Butterfly needle**



- | Blood collection needle, butterfly type.
- | Consist of stainless steel needles (veni-puncturing end and puncturing end with rubber protector), and connecting tube.
- | Universal fitting with plastic needle holders
- | International color-coded
- | Medical grade stainless steel needles
- | Tube length: 12 inches.
- | Gauges: 18G-25G available.
- | EO sterile, Non toxic, non pyrogenic
- | Validity: 3 years

Cat No.	Product	Gauge x Length	Color	Packing info
N-N3M1	Blood collection needle, butterfly type, sterile	20G x 1"	Yellow	50pcs/bag, 1200pcs/ctn
N-N3P1	Blood collection needle, butterfly type, sterile	21G x 3/4"	Green	50pcs/bag, 1200pcs/ctn
N-N3Q1	Blood collection needle, butterfly type, sterile	21G x 1"	Green	50pcs/bag, 1200pcs/ctn
N-N3S1	Blood collection needle, butterfly type, sterile	22G x 3/4"	Black	50pcs/bag, 1200pcs/ctn
N-N3T1	Blood collection needle, butterfly type, sterile	22G x 1"	Black	50pcs/bag, 1200pcs/ctn
N-N3U1	Blood collection needle, butterfly type, sterile	23G x 3/4"	Blue	50pcs/bag, 1200pcs/ctn
N-N3V1	Blood collection needle, butterfly type, sterile	23G x 1"	Blue	50pcs/bag, 1200pcs/ctn

**Mfg.: Nubenco (Shanghai) Co, Ltd**  
**Shanghai, China**  
**Tel.: +86-21-58967745**  
**[www.nubenco.com.cn](http://www.nubenco.com.cn)**