


Fiberglass Orthopedic Casting Tape

Nucast[®]

Fiberglass bandage



- I Fast Setting time 3-5 minutes
 - I Size: 2", 3", 4", 5", 6"
 - I Length: 4 yards
 - I Other sizes available
- I Quickly and easily conform to all shape
- I Low tack resin, allow easy handling during molding
- I Foil bag packed individually, preserving the products from humidity
 - I Box or bag available for double package
 - I Full information printed in the Foil bag
 - I Bar-code available
 - I Non Sterile
 - I Country of Origin: China

More...

- U Complete immobilization line: Fiberglass bandage, POP bandage, Cast padding, Stockinette
- U Size: 2", 3", 4", 5", 6", 8" x 3-4-5 yards
- U Synthetic stockinette(65% Terylene / 35% Cotton), for lower absorption and greater dryness and better air permeability, preferred over 100% cotton
- U Acrylic cast padding (100% Rayon), acts as an air permeable cushion system controlling temperature changes, preferred over other materials.

Cat No.	Product	Graduation	Carton size(cm)	Packing
N-GS2	Orthopedic fiberglass casting tape	2" x 4 yards	50 x 42 x 30	Box x 12 rolls, 96 rolls/carton
N-GS3	Orthopedic fiberglass casting tape	3" x 4 yards	50 x 50 x 30	Box x 12 rolls, 96 rolls/carton
N-GS4	Orthopedic fiberglass casting tape	4" x 4 yards	50 x 32 x 30	Box x 12 rolls, 48 rolls/carton
N-GS5	Orthopedic fiberglass casting tape	5" x 4 yards	50 x 38 x 30	Box x 6 rolls, 48 rolls/carton
N-GS6	Orthopedic fiberglass casting tape	6" x 4 yards	50 x 42 x 30	Box x 6 rolls, 48 rolls/carton
N-GSVB	Orthopedic fiberglass bandage Ultra, Lightweight	2" x 4 yards	70 x 45 x31	Box x 10 rolls, 100 rolls/carton
N-GSWB	Orthopedic fiberglass bandage Ultra, Lightweight	3" x 4 yards	100 x 45 x38	Box x 10 rolls, 100 rolls/carton
N-GSXB	Orthopedic fiberglass bandage Ultra, Lightweight	4" x 4 yards	100 x 45 x38	Box x 10 rolls, 100 rolls/carton
N-GSYB	Orthopedic fiberglass bandage Ultra, Lightweight	5" x 4 yards	100 x 53x38	Box x 10 rolls, 100 rolls/carton
N-GSZB	Orthopedic fiberglass bandage Ultra, Lightweight	6" x 4 yards	82 x 43 x50	Box x 10 rolls, 100 rolls/carton
N-G72	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	2" x 3 yards	63 x 33 x25	Box x 12 rolls, 288 rolls/carton
N-G73	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	3" x 3 yards	48 x 33 x25	Box x 12 rolls, 144 rolls/carton
N-G74	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	4" x 3 yards	63 x 33 x25	Box x 12 rolls, 144 rolls/carton
N-G75	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	5" x 3 yards	63 x 33 x25	Box x 12 rolls, 144 rolls/carton

N-G76	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	6" x 3 yards	63 x 33 x38	Box x 6 rolls, 144 rolls/carton
N-G77	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	8" x 3 yards	63 x 33 x50	Box x 6 rolls, 288 rolls/carton
N-GRW	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	3" x 4 yards		Bag x 12 rolls, 144 rolls/carton
N-GRX	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	4" x 4 yards		Bag x 12 rolls, 144 rolls/carton
N-GRY	Orthopedic Cast Padding, Synthetic, 75g, 50% Rayon, 50% Cotton	5" x 4 yards		Bag x 6 rolls, 144 rolls/carton
N-GTV	Cast padding, acrylic, 100% Rayon	2" x 3 yards	50 x 38 x 24	individual, 192 rolls/carton
N-GTW	Cast padding, acrylic, 100% Rayon	3" x 3 yards	50 x 38 x 34	individual, 192 rolls/carton
N-GTX	Cast padding, acrylic, 100% Rayon	4" x 3 yards	50 x 38 x 24	individual, 96 rolls/carton
N-GTW	Cast padding, acrylic, 100% Rayon	5" x 3 yards	50 x 38 x 24	individual, 96 rolls/carton
N-GTZ	Cast padding, acrylic, 100% Rayon	6" x 3 yards	50 x 38 x 34	individual, 96 rolls/carton
N-G923	Stockinette (tubular bandage), Synthetic, 65% Terylene, 35% Cotton	2" x 25 yards	42 x 42 x33	20 rolls/carton
N-G931	Stockinette (tubular bandage), Synthetic, 65% Terylene, 35% Cotton	3" x 25 yards	42 x 42 x33	16 rolls/carton
N-G924	Stockinette (tubular bandage), Synthetic, 65% Terylene, 35% Cotton	4" x 25 yards	42 x 42 x33	12 rolls/carton
N-G941	Stockinette (tubular bandage), Synthetic, 65% Terylene, 35% Cotton	5" x 25 yards	42 x 42 x33	8 rolls/carton
N-G925	Stockinette (tubular bandage), Synthetic, 65% Terylene, 35% Cotton	6" x 25 yards	42 x 42 x33	8 rolls/carton

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

Fax: +86-21-58967742

[Email:info@nubenco.com.cn](mailto:info@nubenco.com.cn)

www.nubenco.com.cn