

## Surgical gloves/Latex gloves

## Sterile

Extra protection for medical staff and patients in aseptic procedures



- | 100% latex
- | Low protein content, hypoallergenic
- | High resistance and elasticity
- | Textured or smooth
- | High sensitivity
- | Lightly powdered with corn starch
- | Powdered free available
- | Validity: 3 years
- | Packing: individually sealed per pair

Cat. No.	Product	Packing
N-L1A	Sterile surgical glove, latex, smooth, lightly powdered, #6.0	Box of 50 pairs, 500 pairs / carton
N-L1A1	Sterile surgical glove, latex, textured, lightly powdered, #6.0	Box of 50 pairs, 500 pairs / carton
N-L1B	Sterile surgical glove, latex, smooth, lightly powdered, #6.5	Box of 50 pairs, 500 pairs / carton
N-L1B1	Sterile surgical glove, latex, textured, lightly powdered, #6.5	Box of 50 pairs, 500 pairs / carton
N-L1C	Sterile surgical glove, latex, smooth, lightly powdered, #7.0	Box of 50 pairs, 500 pairs / carton
N-L1C1	Sterile surgical glove, latex, textured, lightly powdered, #7.0	Box of 50 pairs, 500 pairs / carton
N-L1D	Sterile surgical glove, latex, smooth, lightly powdered, #7.5	Box of 50 pairs, 500 pairs / carton
N-L1D1	Sterile surgical glove, latex, textured, lightly powdered, #7.5	Box of 50 pairs, 500 pairs / carton
N-L1E	Sterile surgical glove, latex, smooth, lightly powdered, #8.0	Box of 50 pairs, 500 pairs / carton
N-L1E1	Sterile surgical glove, latex, textured, lightly powdered, #8.0	Box of 50 pairs, 500 pairs / carton
N-L1F	Sterile surgical glove, latex, smooth, lightly powdered, #8.5	Box of 50 pairs, 500 pairs / carton
N-L1F1	Sterile surgical glove, latex, textured, lightly powdered, #8.5	Box of 50 pairs, 500 pairs / carton

N-L5A	Sterile surgical glove, latex, non powdered, #6.0	Box of 50 pairs, 500 pairs / carton
N-L5B	Sterile surgical glove, latex, non powdered, #6.5	Box of 50 pairs, 500 pairs / carton
N-L5C	Sterile surgical glove, latex, non powdered, #7.0	Box of 50 pairs, 500 pairs / carton
N-L5D	Sterile surgical glove, latex, non powdered, #7.5	Box of 50 pairs, 500 pairs / carton
N-L5E	Sterile surgical glove, latex, non powdered, #8.0	Box of 50 pairs, 500 pairs / carton
N-L5F	Sterile surgical glove, latex, non powdered, #8.5	Box of 50 pairs, 500 pairs / carton

**Mfg.: Nubenco (Shanghai) Co, Ltd**

**Shanghai, China**

**Tel.: +86-21-58967745**

**Fax: +86-21-58967742**

**[Email:info@nubenco.com.cn](mailto:info@nubenco.com.cn)**

**[www.nubenco.com.cn](http://www.nubenco.com.cn)**