

## Surgical Sutures



- Ø Suture materials: POLYGLYCOLIC ACID , CHROMIC(Catgut), PLAIN(Catgut), SILK(Braided), NYLON(monofilament), POLYPROPYLENE (monofilament), POLYESTER (Braided), STAINLESS STEEL WIRE(monofilament), LIGATURE
- Ø Needle shapes: 1/2 circle, 3/8 circle, 1/4 circle, 5/8 circle
- Ø Needle point type: taper point, blunt taper point, cutting edge, reverse cutting edge, tapercut, micro-point spatula curved

Ø Sterile: EO

Ø Validity: 5 years




POLYGLYCOLIC ACID PG



CHROMIC CATGUT CC



PLAIN CATGUT CP



POLYPROPYLENE PM




POLYESTER BRAIDED PB





SILK BRAIDED SK





















STAINLESS STEEL WIRE SW




NYLON NL

THREAD DIAMETER		THREAD LENGTH		NEEDLE CURVATURE	
10	100	45CM	45		STRAIGHT S1
9/0	90	50CM	50		1/2 CIRCLE C1

8/0	80	60CM	60		1/2 CIRCLE (DOUBLE)	C2
7/0	70	70CM	70		1/4 CIRCLE	C3
6/0	60	75CM	75		1/4 CIRCLE (DOUBLE)	C4
5/0	50	77CM	77		3/8 CIRCLE	C5
4/0	40	90CM	90		3/8 CIRCLE (DOUBLE)	C6
3/0	30	100CM	100		5/8 CIRCLE	C7
2/0	20	125CM	125		1/2 CURVE	C8
1/0	10	150CM	150		LOOP ROUND	L1
0	00					
1	01					
2	02					
3	03					

NEEDLE POINT TYPE	CROSS SECTION	NEEDLE LENGTH
	ROUND BODIED R1	6MM 06
	ROUND BODIED (HEAVY) R2	8MM 08
	CUTTING CV1	10MM 10
	CUTTING (HEAVY) CV2	12MM 12
	REVERSE CUTTING RV1	13MM 13
	REVERSE CUTTING (HEAVY) RV2	15MM 15
	TAPERCUT T1	16MM 16
	TAPERCUT (HEAVY) T2	18MM 18
	BLUNT TAPER POINT T3	20MM 20
	MICRO-POINT SPATULA CURVED M1	30MM 30
	MICRO-POINT SPATULA CURVED M2	35MM 35

	HEAVY	
	TAPER POINT	T4

40MM	40
45MM	45

**Mfg.: Nubenco (Shanghai) Co, Ltd**  
**Shanghai, China**

**Tel.: +86-21-58967745**

**Fax.: +86-21-58967742**

**Email:info@nubenco.com.cn**

**www.nubenco.com.cn**