

Ultra accuracy in dosage and safety

Disposable Syringe



- Ø 3 pcs: barrel, plunger, tip
- Ø With or without needles optional
- Ø Attached high grade stainless steel needle and latex-free optional
- Ø Luer lock or luer slip
- Ø Breakable plunger-NO REUSE
- Ø Siliconized for easy and smooth plunger slide
- Ø Positive plunger stop

- Ø Permanent engraved fine scale marking on barrel
- Ø Grid plunger top for controlled injection
- Ø EO sterile
- Ø individual blister or PE packed available
- Ø Validity: 5 years

CE 0197

ISO13485

Cat No.	Product	Graduation	Gauge x Length	Packing
N-S4H	3ml, luer lock, no needle	3ml		Box x 100pcs, 1800pcs/carton
N-S1HQ	3ml, luer lock, needle 21G X1 1/2"	3ml	21G X1 1/2"	Box x 100pcs, 1800pcs/carton
N-S1HT	3ml, luer lock, needle 22G X1 1/2"	3ml	22G x 1 1/2"	Box x 100pcs, 1800pcs/carton
N-S1HV	5ml, luer lock, needle 23G X1"	3ml	23G X1"	Box x 100pcs, 1800pcs/carton
N-S4I	5ml, luer lock, no needle	5ml		Box x 100pcs, 1200pcs/carton
N-S1IW	5ml, luer lock, needle 23G X1 1/2"	5ml	23G X1 1/2"	Box x 100pcs, 1200pcs/carton
N-S1IQ	5ml, luer lock, needle 21G X1 1/2"	5ml	21G X1 1/2"	Box x 100pcs, 1200pcs/carton

N-S4J	10ml, luer lock, no needle	10ml		Box x 100pcs, 1200pcs/carton
N-S1JQ	10ml, luer lock, needle 21G X1 1/2"	10ml	21G x 1 1/2"	Box x 100pcs, 1200pcs/carton
N-S1KQ	20ml, luer lock, needle 21G X1 1/2"	20ml	21G x 1 1/2"	Box x 50pcs, 600pcs/carton
N-SA4K	20ml, luer lock, no needle	20ml		Box x 50pcs, 600pcs/carton
N-S3M	30ml, luer lock, no needle	30ml		Box x 50pcs, 400pcs/carton
N-S3P	50ml, luer lock, no needle	50ml		Box x 50pcs, 400pcs/carton
N-S4P	50ml, luer slip, no needle	50ml		Box x 50pcs, 400pcs/carton
N-S3R	60ml, luer lock, no needle	60ml		Box x 50pcs, 400pcs/carton
N-S4R	60ml, luer slip, no needle	60ml		Box x 50pcs, 400pcs/carton

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

Fax.: +86-21-58967742

Email:info@nubenco.com.cn

www.nubenco.com.cn