

Sterilization Paper
SteriPaP[®]
Medical Grade Paper



- I Applicable for vapor and ethylene oxide gas.**
- I One side medical grade paper and other side transparent film Polyester-Polyethylene.**
- I Flat or gusseted rolls and pouches in all size available.**

Easy to use, clear view of products inside

Longer sterility after sterilization

60g/m² good density paper

Indicators along edges of paper roll

Repels liquid

Antimicrobial barrier

Seals well

No ETO residue

Complies with EN868

| Cat. No | Size (metric system) | Roll type | Packaging |
|---------|-------------------------|----------------|---------------|
| N-PF1A | 50mm x 200m | FLAT (ROLL) | 6 rolls / box |
| N-PF1B | 75mm x 200m | FLAT (ROLL) | 4 rolls / box |
| N-PF1C | 100mm x 200m | FLAT (ROLL) | 3 rolls / box |
| N-PF1D | 150mm x 200m | FLAT (ROLL) | 2 rolls / box |
| N-PF1E | 200mm x 200m | FLAT (ROLL) | 1 roll / box |
| N-PF1F | 250mm x 200m | FLAT (ROLL) | 1 roll / box |
| N-PF1G | 300mm x 200m | FLAT (ROLL) | 1 roll / box |
| N-PF1H | 350mm x 200m | FLAT (ROLL) | 1 roll / box |
| N-PF1I | 400mm x 200m | FLAT (ROLL) | 1 roll / box |
| N-PG1C | 100mm x 100m | GUSSETED(ROLL) | 2 rolls / box |
| N-PG1D | 150mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PG1E | 200mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PG1F | 250mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PG1G | 300mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PG1H | 350mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PG1I | 400mm x 100m | GUSSETED(ROLL) | 1 roll / box |
| N-PSP | 40 x 50cm | FLAT (POUCH) | 1 each |
| N-PSY | 10 x 20cm | FLAT (POUCH) | 1 each |
| N-PSZ | 20 x 35cm | FLAT (POUCH) | 1 each |

Mfg.: Nubenco (Shanghai) Co, Ltd

Shanghai, China

Tel.: +86-21-58967745

Fax: +86-21-58967742

[Email:info@nubenco.com.cn](mailto:info@nubenco.com.cn)

www.nubenco.com.cn